# **Sheet B**


# **INSTRUCTIONS FOR WEARING AND ADJUSTMENTS**

NOTE: Keep in mind that the Multi-Mask is a personal product so you will get used to all the parts in just a few days.

### **Tightening the STRAPS**



Raise this part slightly and pull To tighten the mask, pull the the strap as required. Don't raise too much or it can break.



strap backwards.



To loosen the mask, pull the strap forward.

Count the lines on the strap and make sure both sides are equal.

You can even cut off the extra lengths of the strap if you want. But please do this after a few weeks of use so that you are sure, and still leave about 3 cm extra to be on the safe side.

#### Putting on the Multi-Mask













For everyday wearing, just use the rear elastic band to loosen for wearing or removing, and for tightening the mask on the face. All other straps don't need to be changed ever again.

## **Tightening the ELASTIC BAND of the Headband**



The elastic strap sits in the smaller hole and remains locked with the beads.



Simply pull the beaded silicone cord towards the bigger hole.



it by moving it in this bigger hole.



You can now loosen or tighten At desired tightness, simply lock it by pulling it into the smaller hole.

#### Putting on the Shield on the Multi-Mask







The shield slides onto the main cup. You will notice a groove around the opening on the body of the main cup, and the shield slides into the groove. This will take a little practice to slide easily.

Make sure the small notch on the inner cut out of the shield sits on the tiny projection on the main cup, this ensures it will not slip sideways. To remove the shield, hold the mask and just slide it off upwards.

### **Tightening the top of the HEADBAND**



part on the top of the headband.



Note this small locking plastic For tightening the headband, simply push the 2 sides towards each other.



The plastic tab is designed to automatically slip lock when it is being tightened.

## **Loosening the top of the HEADBAND**



With index finger hold the end Once the lock is released, use projection and with thumb pull your other hand to pull the the small tab up to release the other side of the headband lock.



apart.



This will make both the parts move apart, loosening the headband.

Tightening the headband is easy but to loosen, it is deliberately made difficult to avoid it from accidentally slipping.